

Standard Operating Procedure:  
**Serum & Plasma Sample Clarification Via Centrifugation**

**1. Distribution:**

<b>Personnel</b>	<b>Training Date</b>	<b>SOP Version</b>	<b>Signature</b>

Standard Operating Procedure:  
**Serum & Plasma Sample Clarification Via Centrifugation**

Date:	02-21-2020		Effective Date:  02-21-2020
Written By:	Danielle Parent, MLS		
Approvals:			
	Laboratory Director	Laboratory Coordinator	
Related Documents & attachments:			
Revision History:			
<u>Rev #</u>	<u>Effective Date</u>	<u>Description of Change</u>	

**2. Background and Principle:**

- 2.1 Background:** In accordance to the Kit Instructions for the Simoa N4PA Advantage Kit (Catalog Number: 102153), excess fibrin, lipids and particulate matter may interfere with the accuracy of the assay. The following procedure efficiently removes these interfering substances.
- 2.2 Principle:** High speed centrifugation will separate lipids (top layer) and fibrin/particulate matter (pellet formation) such that a clarified serum or plasma sample can be plated for use on the Quanterix HD-1 or HD-X analyzer.

**3. Specimen Handling and Collection:**

**3.1 Specimen Handling:**

- 3.1.1** Blood-borne pathogen safety training (offered by UVM ESF) must be completed before any work in the laboratory commences. **Use Universal Precautions and treat blood and blood products as potentially infectious materials. It is not known if specimens contain HBV, HIV, and other blood borne pathogens.**
- 3.1.2** Wear proper protective equipment, lab coat and gloves.
- 3.1.3** Dispose all tips and materials that come in contact with biological agents into proper biological waste containers.

Standard Operating Procedure:  
**Serum & Plasma Sample Clarification Via Centrifugation**

**3.2 Specimen Collection:**

**3.2.1 Serum:** Use a serum separator tube (SST) and allow samples to clot for 30 minutes before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at < -20 o C. Avoid repeated freeze-thaw cycles.

**3.2.2 Plasma:** Collect plasma using EDTA or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at < -20 o C. Avoid repeated freeze-thaw cycles.

**4. Reagents, Recipes and Equipment**

**4.1 Reagents:** 10% Bleach for work area decontamination

**4.2 Equipment:**

**4.2.1 Vacuum Apparatus**

- Bench Top Built in or equivalent

**4.2.2 1L - 2L Vacuum flask**

- Thermo Fisher Catalog DS4101-1000 or equivalent

**4.2.3 One hole Stopper to fit vacuum flask**

- Fisher Scientific Catalog 14-135M (Stopper size 8)

**4.2.4 Tygon Tubing**

- ¼ - ⅜ inch or equivalent to fit Vacuum Flask.
- Thermo Fisher Catalog 8702-0065 or equivalent

**4.2.5 Vacuum VacTrap (Alternate to individual pieces listed above)**

- Fisher Scientific Catalog 50-148-9320

**4.2.6 Adjustable Pipette (50 – 300uL & 100 – 1000uL)**

**4.2.7 200 – 300uL Pipette Tips**

**4.2.8 1000uL Pipette Tips**

**4.2.9 Vortexer/Mixer**

- Thermo Fisher Catalog 88880017TS or equivalent.

Standard Operating Procedure:  
**Serum & Plasma Sample Clarification Via Centrifugation**

**4.2.10** Centrifuge/Micro-centrifuge:

- Thermo Fisher Catalog 75002437 or equivalent, capable of 14,000g.

**4.2.11** Bright Work Lamp

- Thermo Fisher Catalog 11-990-108 or equivalent

**4.2.12** Cryovial, 12x75 Test Tube or 15mL Conical or other suitable temporary container to mix clarified sample.

- 1.5mL Conical Tube
  - 4.2.12..1** Fisher Scientific Catalog 50-754-1456
- 12x75 Polystyrene Test Tubes
  - 4.2.12..1** Fisher Scientific Catalog 22-171-604
- 15mL Conical Centrifuge Tube
  - 4.2.12..1** Fisher Scientific Catalog 14-959-53A

**5. Procedural Steps:**

**5.1.1** Method 1: Low Sample Volume

- **Step 1:** Thaw sample at 37°C for 5-7 minutes
- **Step 2:** Vortex sample to mix well
- **Step 3:** Centrifuge sample(s) at 14,000g for 10 minutes
- **Step 4:** Carefully remove the vial from the centrifuge such that the lipid layer and pellet formed is not disturbed.
- **Step 5:** Carefully insert pipette into the clarified serum or plasma. Take care to minimally disturb any lipid layer and/or pellet. Tilting the cryo may be helpful in exposing clarified serum.
- **Step 6:** Aliquot to desired receptacle: cryovial, HD-1/HD-X plate.
- **Step 7:** Freeze aliquots at -80°C until further use. And/or Run on HD-1/HD-X analyzer

Standard Operating Procedure:  
**Serum & Plasma Sample Clarification Via Centrifugation**

**5.1.2 Method 2:** High volume and/or multiple aliquot procedure

- **Step 1:** Thaw sample at 37°C for 5-7 minutes, or until sample is fully thawed.
- **Step 2:** Vortex Sample to mix well.
- **Step 3:** Centrifuge sample(s) at 14,000g for 10 minutes
- **Step 4:** Place 200uL or 300uL tip onto the end of the vacuum tubing
- **Step 5:** Slowly and carefully insert vacuum tip into the centrifuged sample and remove lipid layer. (Expect to lose ~200uL – 250uL depending on amount of lipids). Replace tip between samples.
- **Step 6:** Carefully remove all clarified serum or plasma to a NEW temporary vial (cryo-vial, 12x75 tube, 15mL conical tube etc.)
- **Step 7:** Ensure that serum or plasma is visually clear of lipids or other debris or particulates. If any debris remains, centrifuge again at 14,000g for 10 minutes.
- **Step 8:** Vortex clarified serum or plasma
- **Step 9:** Transfer clarified serum or plasma into desired daughter aliquots.
- **Step 10:** Cap and freeze aliquots at -80°C until further use.

Standard Operating Procedure:  
**Serum & Plasma Sample Clarification Via Centrifugation**

**Certificate of Review**

<b>Title:</b>			
<b>Prepared By:</b>		<b>Date Prepared:</b>	_/_/____
<b>Distribution:</b>		<b>Location:</b>	
<b>Supersedes:</b>		<b>Date Initiated:</b>	_/_/____
<b>Comments:</b>		<b>By:</b>	

<b>Date Reviewed:</b>	_/_/____	<b>By:</b>	
<b>Comments:</b>			

<b>Date Reviewed:</b>	_/_/____	<b>By:</b>	
<b>Comments:</b>			

<b>Date Reviewed:</b>	_/_/____	<b>By:</b>	
<b>Comments:</b>			

<b>Date Reviewed:</b>	_/_/____	<b>By:</b>	
<b>Comments:</b>			

<b>Date Reviewed:</b>	_/_/____	<b>By:</b>	
<b>Comments:</b>			

<b>Date Reviewed:</b>	_/_/____	<b>By:</b>	
<b>Comments:</b>			